

NEC YORK

Kit 3 sizes one-piece cervical collars



Nec York cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features


- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- 3 sizes (paediatric, neckless, regular)
- Storage bag

Standard equipment

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Technical data

Sizes	2 (Paediatric) 3 (Neckless) 5 (Regular)
Materials	PE, PA
Kit weight (with bag)	430 ± 50 g
Dimension A (± 5mm)	125 mm (Paediatric) 135 mm (Neckless) 145 mm (Regular)
Dimension B (±5 mm)	445 mm (Paediatric) 555 mm (Neckless) 555 mm (Regular)
Dimension C (±5 mm)	100 mm (Paediatric) 95 mm (Neckless) 115 mm (Regular)



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

QC19013C

NEC YORK KIT 3SIZES(2,3,5)WITH BAG1PIECE C.COLL.

SET consisting of the following MD compliant with UE Reg. 2017/745

QC19002C – NEC YORK SZ.2 - 1 PIECE RIGID CERVICAL COLLAR PEDI	Reg. n. 83084	CND Y060312
QC19020C – NEC YORK SIZE 3 (NECKLESS) - 1PIECE CERV. COLLAR	Reg. n. 83086	CND Y060312
QC19022C – NEC YORK SIZE 5 (REGULAR) - 1PIECE CERV. COLLAR	Reg. n. 83090	CND Y060312

Rev.0 (31/03/2021)

 UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it